

NED University of Engineering & Technology, Karachi

Department of Software Engineering



Formal Methods in Software Engineering(SE-313)
ASSIGMENT

Professor: Dr.Mustafa Latif

Title: Ambulatory Blood Pressure Monitor System

Group: B7

Members:

Khawar Khan(SE-21093)

Muhammad Zubair(SE-21094)

DATE: Jan 11, 2024

INDEX

S.No	Section	Page Number
01	Scope	Page#1
02	4+1 view model	Page#2
03	VDM Specifications and Java Code	Page#7
04	VDM specification of ABPM	Page#7
05	Java Code	Page#13
06	Testing class test cases summary	Page#27
07	Tester class Java Code	Page#31

Ambulatory Blood Pressure Monitoring System

I. SCOPE:

The project focuses on developing an Ambulatory Blood Pressure Monitoring (ABPM) system designed to provide continuous and accurate monitoring of blood pressure levels. The primary goal is to assist healthcare professionals in diagnosing and managing various blood pressure conditions, including Normal, Elevated, Hypertension Stage 1, Hypertension Stage 2, and Hypertensive Crisis.

Here are the general categories for blood pressure levels:

- Low Blood Pressure: Systolic < 90 mm Hg and Diastolic < 60 mm Hg
- Normal: Systolic < 120 mm Hg and Diastolic < 80 mm Hg
- Elevated: Systolic 120-129 mm Hg and Diastolic < 80 mm Hg
- Hypertension Stage 1: Systolic 130-139 mm Hg or Diastolic 80-89 mm Hg
- Hypertension Stage 2: Systolic \geq 140 mm Hg or Diastolic \geq 90 mm Hg
- Hypertensive Crisis: Systolic > 180 mm Hg and/or Diastolic > 120 mm Hg

The system's core functionality lies in its ability to perform automatic and regular blood pressure measurements at predetermined intervals throughout the regular period. These measurements will be categorized based on internationally recognized standards, allowing for the classification of blood pressure readings into distinct categories.

The integrity of the ABPM system is paramount, representing a cornerstone of trust for healthcare professionals in making critical decisions for patient care. The system's high level of integrity is manifested in its unwavering commitment to accurate and continuous blood pressure monitoring. Through robust algorithms and adherence to international standards, the system ensures precise categorization of blood pressure readings, contributing to an enhanced level of diagnostic precision. This reliability becomes the linchpin for healthcare practitioners, empowering them with trustworthy data that directly influences treatment plans and interventions. In this way, the system's commitment to integrity not only elevates diagnostic accuracy but also plays a central role in upholding patient safety throughout the monitoring process.


The anticipated outcomes of this ABPM system include improved diagnostic accuracy, personalized treatment plans, and continuous patient care. By providing

Ambulatory Blood Pressure Monitoring System

healthcare professionals with a robust tool for monitoring blood pressure, the system aims to contribute to better health outcomes and increased patient engagement in the management of their blood pressure conditions.

II. 4+1 VIEW MODEL


1. Logical View



Activate Windows


Ambulatory Blood Pressure Monitoring System

2. Process View



Ambulatory Blood Pressure Monitoring System


3. Physical View Deployment Diagram



Ambulatory Blood Pressure Monitoring System


4. Development View

a. Package Diagram




Ambulatory Blood Pressure Monitoring System

b. Component Diagram



Ambulatory Blood Pressure Monitoring System

5. Scenarios View



III. VDM-SL SPECIFICATION OF THE SYSTEM & JAVA CODE

IV. VDM-Specification of APBM

types

SIGNAL = <SILENT> | <VIBRATE> | <VIBRATE_AND_ALERT>

BLOODPRESSURELEVEL = <LOW_BLOOD_PRESSURE> | <NORMAL> |
<ELEVATED> | <HBP_STAGE_1> | <HBP_STAGE_2> |
<HYPERTENSIVE_CRISIS>

Values

MIN_BLOOD_PRESSURE_LEVEL: Z = 1

MAX_BLOOD_PRESSURE_LEVEL: Z = 400

State *AmbulatoryBloodPressureController* of

systolic: [Z]

diastolic: [Z]

currentBloodPressureLevel : [BloodPressureLevel]

Ambulatory Blood Pressure Monitoring System

currentSignal : [Signal]
firstname: [seq of char]
lastname: [seq of char]
age:[Z]

inv mk-AmbulatoryBloodPressureController (s,d) Δ ($MIN_BLOOD_PRESSURE_LEVEL \leq s \leq MAX_BLOOD_PRESSURE_LEVEL \vee s=0$) \wedge ($MIN_BLOOD_PRESSURE_LEVEL \leq d \leq MAX_BLOOD_PRESSURE_LEVEL \vee d=0$)

init mk-AmbulatoryBloodPressureController (s,d,fn,ln,a) Δ ($s=0$) \wedge ($d=0$) \wedge ($fn=|$) \wedge ($an=|$) \wedge ($ln=|$)

end

operations

inRange(val: Z) result: B

pre TRUE

post result $\Leftrightarrow MIN_BLOOD_PRESSURE_LEVEL \leq val \leq MAX_BLOOD_PRESSURE_LEVEL$

setInitialSystolicPressure(pressureIn:Z)

ext wr systolic: [Z]

pre inRange(pressureIn) \wedge systolic=0

post systolic = pressureIn

setInitialDiastolicPressure(pressureIn: Z)

ext wr diastolic: [Z]

pre inRange(pressureIn) \wedge diastolic=0

post diastolic = pressureIn

getCurrentBloodPressureLevel() : bloodPressureLevelOut : BloodPressureLevel

ext rd currentBloodPressureLevel

pre currentBloodPressureLevel \neq nil

post bloodPressureLevelOut = currentBloodPressureLevel

getCurrentSignal() : signalOut : Signal

ext rd currentSignal

pre currentSignal \neq nil

post signalOut = currentSignal

changeSystolicPressure(pressureIn) bloodPressureLevelOut: BloodPressureLevel

ext wr systolic: [Z]

pre inRange(pressureIn)

post systolic = pressureIn \wedge bloodPressureLevelOut = currentBloodPressure

Ambulatory Blood Pressure Monitoring System

changeDiastolicPressure(pressureIn) bloodPressureLevelOut: BloodPressureLevel

ext wr diastolic: [Z]

pre inRange(pressureIn)

post diastolic = pressureIn \wedge bloodPressureLevelOut = currentBloodPressure

determineBloodPressureLevel() currentBloodPressureLevelOut : BloodPressureLevel

ext rd systolic: [Z]

rd diastolic: [Z]

ext wr currentBloodPressureLevel

pre inRange(systolic) \wedge systolic $\neq 0$ \wedge inRange(diastolic) \wedge diastolic $\neq 0$

post ((systolic < 90 \wedge diastolic < 60 \wedge currentBloodPressureLevel = <LOW_BLOOD_PRESSURE> \wedge currentBloodPressureLevelOut = <LOW_BLOOD_PRESSURE>) \vee (90 \leq systolic \leq 120 \wedge 60 \leq diastolic \leq 80 \wedge currentBloodPressureLevel = <NORMAL> \wedge currentBloodPressureLevelOut = <NORMAL>) \vee (120 \leq systolic \leq 130 \wedge diastolic \leq 80 \wedge currentBloodPressureLevel = <ELEVATED> \wedge currentBloodPressureLevelOut = <ELEVATED>) \vee (130 < systolic \leq 140 \wedge 80 \leq diastolic \leq 90 \wedge currentBloodPressureLevel = <HBP_STAGE_1> \wedge currentBloodPressureLevelOut = <HBP_STAGE_1>) \vee (systolic > 140 \wedge diastolic > 90 \wedge currentBloodPressureLevel = <HYPERTENSIVE_CRISIS> \wedge currentBloodPressureLevelOut = <HYPERTENSIVE_CRISIS>) \vee (currentBloodPressureLevel = <HBP_STAGE_2> \wedge currentBloodPressureLevelOut = <HBP_STAGE_2>))

determineSignal() signalOut : Signal

ext rd systolic: [Z]

rd diastolic: [Z]

ext wr currentSignal: [Signal]

pre inRange(systolic) \wedge systolic $\neq 0$ \wedge inRange(diastolic) \wedge diastolic $\neq 0$

post ((systolic < 90 \wedge diastolic < 60 \wedge currentSignal = <LOW_BLOOD_PRESSURE> \wedge signalOut = <ALERT_AND_VIBRATE>) \vee (90 \leq systolic \leq 120 \wedge 60 \leq diastolic \leq 80 \wedge currentSignal = <SILENT> \wedge signalOut = <SILENT>) \vee (120 \leq systolic \leq 130 \wedge diastolic \leq 80 \wedge currentSignal = <VIBRATE> \wedge signalOut = <VIBRATE>) \vee (130 < systolic \leq 140 \wedge 80 \leq diastolic \leq 90 \wedge currentSignal = <ALERT_AND_VIBRATE> \wedge signalOut = <ALERT_AND_VIBRATE>) \wedge signalOut = <VIBRATE_AND_ALERT>) \vee (systolic > 140 \wedge diastolic > 90 \wedge currentSignal = <VIBRATE_AND_ALERT> \wedge signalOut = <VIBRATE_AND_ALERT>) \vee (currentSignal = <VIBRATE_AND_ALERT> \wedge signalOut = <VIBRATE_AND_ALERT>))

getSystolicPressure() pressureOut: [Z]

ext rd systolic: [Z]

pre systolic $\neq 0$

Ambulatory Blood Pressure Monitoring System

post pressureOut = systolic

getDiastolicPressure() pressureOut: [Z]

ext rd diastolic: [Z]

pre diastolic \neq 0

post pressureOut = diastolic

VDM Specification	Java Code
<p>types</p> <p><i>SIGNAL</i> = <SILENT> <VIBRATE> <VIBRATE_AND_ALERT> <i>BLOODPRESSURELEVEL</i> = <LOW_BLOOD_PRESSURE> <NORMAL> <ELEVATED> <HBP_STAGE_1> <HBP_STAGE_2> <HYPERTENSIVE_CRISIS></p>	<pre>enum BloodPressureLevel { LOW_BLOOD_PRESSURE, NORMAL, ELEVATED, HBP_STAGE_1, HBP_STAGE_2, HYPERTENSIVE_CRISIS } enum SIGNAL { SILENT, VIBRATE, VIBRATE_AND_ALERT }</pre>
<p>values</p> <p><i>MIN_BLOOD_PRESSURE_LEVEL</i>: Z = 1 <i>MAX_BLOOD_PRESSURE_LEVEL</i>: Z = 400</p>	<pre>public final int MIN_BLOOD_PRESSURE_LEVEL = 1; public static final int MAX_BLOOD_PRESSURE_LEVEL = 400;</pre>
<p>inv mk-AmbulatoryBloodPressureController (s,d) $\underline{\Delta}$ (<i>MIN_BLOOD_PRESSURE_LEVEL</i> \leq s \leq <i>MAX_BLOOD_PRESSURE_LEVEL</i>) \vee s=0) \wedge (<i>MIN_BLOOD_PRESSURE_LEVEL</i> \leq d \leq</p>	<pre>public boolean inv() { return (((this.systolic >= MIN_BLOOD_PRESSURE_LEVEL &&</pre>

Ambulatory Blood Pressure Monitoring System

$MAX_BLOOD_PRESSURE_LEVEL \vee d=0$	<pre>this.systolic<=MAX_BLOOD_PRESSURE_LEVEL) this.systolic==0) && (((this.diastolic>=MIN_BLOOD_PRESSURE_LEVEL && this.diastolic<=MAX_BLOOD_PRESSURE_LEVEL) this.systolic==0); }</pre>
<hr/> <p>init mk-AmbulatoryBloodPressureController (s,d,fn,ln,a) Δ (s=0) \wedge (d=0) \wedge (fn=) \wedge (an=) \wedge (ln=)</p>	<hr/> <pre>public AmbulatoryBloodPressureController(String fname, String lname, int age) { this.firstname = fname; this.lastname = lname; this.age = age; this.systolic=0; this.diastolic=0; }</pre>
<hr/> <p>inRange(val: Z) result: B pre TRUE post result \Leftrightarrow $MIN_BLOOD_PRESSURE_LEVEL \leq val \leq$ $MAX_BLOOD_PRESSURE_LEVEL$</p>	<hr/> <pre>private boolean inRange(int val) { return val >= MIN_BLOOD_PRESSURE_LEVEL && val <= MAX_BLOOD_PRESSURE_LEVEL; }</pre>
<hr/> <p>setInitialSystolicPressure(pressureIn:Z) ext wr systolic: [Z] pre inRange(pressureIn) \wedge systolic=0 post systolic = pressureIn</p>	<hr/> <pre>public void setInitialSystolicPressure(int pressureIn) { VDM.preTest(inRange(pressureIn)&& systolic==0); this.systolic = pressureIn; VDM.invTest(this); }</pre>
<hr/> <p>setInitialDiastolicPressure(pressureIn: Z) ext wr diastolic: [Z] pre inRange(pressureIn) \wedge diastolic=0 post diastolic = pressureIn</p>	<hr/> <pre>public void setInitialDiastolicPressure(int pressureIn) { VDM.preTest(inRange(pressureIn)&& diastolic==0); this.diastolic = pressureIn; VDM.invTest(this); }</pre>
<hr/> <p>determineBloodPressureLevel() currentBloodPressureLevelOut :</p>	<hr/> <pre>public BloodPressureLevel</pre>

Ambulatory Blood Pressure Monitoring System

```

BloodPressureLevel
ext rd systolic: [Z]
    rd diastolic: [Z]
ext wr currentBloodPressureLevel
pre    inRange(systolic)  $\wedge$  systolic  $\neq$  0  $\wedge$ 
inRange(diastolic)  $\wedge$  diastolic  $\neq$  0
post  ((systolic < 90  $\wedge$  diastolic < 60  $\wedge$ 
currentBloodPressureLevel =
    <LOW_BLOOD_PRESSURE>  $\wedge$ 
currentBloodPressureLevelOut =
    <LOW_BLOOD_PRESSURE>)  $\vee$  ( 90  $\leq$ 
systolic  $\leq$  120  $\wedge$  60  $\leq$  diastolic  $\leq$  80  $\wedge$ 
currentBloodPressureLevel =
    <NORMAL>  $\wedge$ 
currentBloodPressureLevelOut =
    <NORMAL>)  $\vee$  ( 120  $\leq$  systolic  $\leq$  130  $\wedge$ 
diastolic  $\leq$  80  $\wedge$ 
currentBloodPressureLevel =
    <ELEVATED>  $\wedge$ 
currentBloodPressureLevelOut =
    <ELEVATED>)  $\vee$  ( 130 < systolic  $\leq$  140
 $\wedge$  80  $\leq$  diastolic  $\leq$  90  $\wedge$ 
currentBloodPressureLevel =
    <HBP_STAGE_1>  $\wedge$ 
currentBloodPressureLevelOut =
    <HBP_STAGE_1>)  $\vee$  ( systolic > 140
 $\wedge$  diastolic > 90  $\wedge$ 
currentBloodPressureLevel =
    <HYPERTENSIVE_CRISIS>  $\wedge$ 
currentBloodPressureLevelOut =
    <HYPERTENSIVE_CRISIS>)  $\vee$ 
(currentBloodPressureLevel =
    <HBP_STAGE_2>  $\wedge$ 
currentBloodPressureLevelOut =
    <HBP_STAGE_2>))

```

```

determineSignal() signalOut : Signal
ext rd systolic: [Z]
    rd diastolic: [Z]
ext wr currentSignal: [Signal]
pre    inRange(systolic)  $\wedge$  systolic  $\neq$  0  $\wedge$ 
inRange(diastolic)  $\wedge$  diastolic  $\neq$  0
post  ((systolic < 90  $\wedge$  diastolic < 60  $\wedge$ 
currentSignal=

```

```

determineBloodPressureLevel() {

    VDM.preTest(inRange(systolic) &&
inRange(diastolic) && systolic!=0 && diastolic!=0);

    if (systolic < 90 && diastolic < 60) {
        currentBloodPressureLevel =
BloodPressureLevel.LOW_BLOOD_PRESSURE;
    } else if (systolic  $\geq$  90 && systolic  $\leq$  120 &&
diastolic  $\geq$  60 && diastolic  $\leq$  80) {
        currentBloodPressureLevel =
BloodPressureLevel.NORMAL;
    } else if (systolic > 120 && systolic  $\leq$  130 &&
diastolic  $\leq$  80) {
        currentBloodPressureLevel =
BloodPressureLevel.ELEVATED;
    } else if (systolic > 130 && systolic  $\leq$  140 &&
diastolic > 80 && diastolic  $\leq$  90) {
        currentBloodPressureLevel =
BloodPressureLevel.HBP_STAGE_1;
    } else if (systolic > 140 && diastolic > 90) {
        currentBloodPressureLevel =
BloodPressureLevel.HYPERTENSIVE_CRISIS;
    } else {
        currentBloodPressureLevel =
BloodPressureLevel.HBP_STAGE_2;
    }

    VDM.invTest(this);

    return currentBloodPressureLevel;
}

```

```

public SIGNAL determineSignal() {
    boolean validInputRange = inRange(systolic) &&
inRange(diastolic) && systolic != 0 && diastolic !=
0;

```

Ambulatory Blood Pressure Monitoring System

<pre> <LOW_BLOOD_PRESSURE> ∧ signalOut = <ALERT_AND_VIBRATE>) v (90 ≤ systolic ≤ 120 ∧ 60 ≤ diastolic ≤ 80 ∧ currentSignal= <SILENT> ∧ signalOut = <SILENT>) v (120 ≤ systolic ≤ 130 ∧ diastolic ≤ 80 ∧ currentSignal= <VIBRATE> ∧ signalOut = <VIBRATE>) v (130 < systolic ≤ 140 ∧ 80 ≤ diastolic ≤ 90 ∧ currentSignal= <ALERT_AND_VIBRATE> ∧ signalOut = <ALERT_AND_VIBRATE>) ∧ signalOut = <VIBRATE_AND_ALERT>) v (systolic > 140 ∧ diastolic > 90 ∧ currentSignal= <VIBRATE_AND_ALERT> ∧ signalOut = <VIBRATE_AND_ALERT>) v (currentSignal = <VIBRATE_AND_ALERT> ∧ signalOut = <VIBRATE_AND_ALERT>)) </pre>	<pre> VDM.preTest(validInputRange); if (systolic < 90 && diastolic < 60) { currentSignal = SIGNAL.VIBRATE_AND_ALERT; } else if (systolic ≥ 90 && systolic ≤ 120 && diastolic ≥ 60 && diastolic ≤ 80) { currentSignal = SIGNAL.SILENT; } else if (systolic > 120 && systolic ≤ 130 && diastolic ≤ 80) { currentSignal = SIGNAL.VIBRATE; } else if (systolic > 130 && systolic ≤ 140 && diastolic > 80 && diastolic ≤ 90) { currentSignal = SIGNAL.VIBRATE_AND_ALERT; } else if (systolic > 140 && diastolic > 90) { currentSignal = SIGNAL.VIBRATE_AND_ALERT; } else { currentSignal = SIGNAL.VIBRATE_AND_ALERT; } VDM.invTest(this); return currentSignal; } </pre>
--	---

V. JAVA CODE OF ABPM system

```

package ABPM;

import FMS.VDM;

enum BloodPressureLevel {

    LOW_BLOOD_PRESSURE,

    NORMAL,

```

Ambulatory Blood Pressure Monitoring System

```
ELEVATED,  
  
HBP_STAGE_1,  
  
HBP_STAGE_2,  
  
HYPERTENSIVE_CRISIS  
}  
  
enum SIGNAL { SILENT, VIBRATE, VIBRATE_AND_ALERT }  
  
interface invChecker {  
  
    public boolean inv();  
  
}  
  
public class AmbulatoryBloodPressureController implements invChecker{  
  
    public String firstname, lastname;  
  
    public int age;  
  
    private int systolic, diastolic;  
  
    private BloodPressureLevel currentBloodPressureLevel = null;  
  
    private SIGNAL currentSignal = null;  
  
  
  
    public final int MIN_BLOOD_PRESSURE_LEVEL = 1;  
  
    public static final int MAX_BLOOD_PRESSURE_LEVEL = 400;  
  
  
  
    public AmbulatoryBloodPressureController(  
  
        String fname, String lname, int age) {  
  
        this.firstname = fname;  
  
        this.lastname = lname;  
  
        this.age = age;  
  
        this.systolic=0;  
  
        this.diastolic=0;
```


Ambulatory Blood Pressure Monitoring System

```
}

@Override

public boolean inv() {

    return ((this.systolic>=MIN_BLOOD_PRESSURE_LEVEL &&
this.systolic<=MAX_BLOOD_PRESSURE_LEVEL)||
this.systolic==0) &&
((this.diastolic>=MIN_BLOOD_PRESSURE_LEVEL && this.diastolic<=MAX_BLOOD_PRESSURE_LEVEL)||
this.systolic==0);

}

private boolean inRange(int val) {

    return val >= MIN_BLOOD_PRESSURE_LEVEL && val <= MAX_BLOOD_PRESSURE_LEVEL;

}

public void setInitialSystolicPressure(int pressureIn) {

    try {

        VDM.preTest(inRange(pressureIn)&& systolic==0);

        this.systolic = pressureIn;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: "

            + "Initial systolic Value given is out of Range");

        this.systolic = 0;

    }

    try{

        VDM.invTest(this);

    } catch(Exception e) {

        System.err.println("Postcondition test failed: "

            + e.getMessage());

    }

}
```

Ambulatory Blood Pressure Monitoring System

```
public void setInitialDiastolicPressure(int pressureIn) {

    try {

        VDM.preTest(inRange(pressureIn));

        this.diastolic = pressureIn;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: "

            + "Initial diastolic Value given is out of Range");

        this.diastolic = 0;

    }

    try{

        VDM.invTest(this);

    } catch(Exception e) {

        System.err.println("Postcondition test failed: "

            + e.getMessage());

    }

}

public BloodPressureLevel changeSystolicPressure(int pressureIn) {

    try {

        VDM.preTest(inRange(pressureIn));

        this.systolic = pressureIn;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: "

            + "Given systolic Value given is out of Range");

        this.systolic = 0;

    }

    try{
```

Ambulatory Blood Pressure Monitoring System

```
VDM.invTest(this);

} catch(Exception e) {

    System.err.println("Postcondition test failed: "

        + e.getMessage());

}

return determineBloodPressureLevel();

}

public BloodPressureLevel changeDiastolicPressure(int pressureIn) {

    try {

        VDM.preTest(inRange(pressureIn));

        this.diastolic = pressureIn;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: "

            + "Given diastolic Value given is out of Range");

        this.diastolic=0;

    }

    try{

        VDM.invTest(this);

    } catch(Exception e) {

        System.err.println("Postcondition test failed: "

            + e.getMessage());

    }

    return determineBloodPressureLevel();

}

public BloodPressureLevel getCurrentBloodPressureLevel() {

    try {
```

Ambulatory Blood Pressure Monitoring System

```
VDM.preTest(currentBloodPressureLevel != null);

return currentBloodPressureLevel;

} catch (VDM.VDMPreconditionException e) {

    System.err.println("Precondition test failed: "

        + "CurrentBloodPressue Level not Intialized");

}

try{

    VDM.invTest(this);

} catch(Exception e) {

    System.err.println("Postcondition test failed: "

        + e.getMessage());

}

return null;

}

public SIGNAL getCurrentSignal() {

    try {

        VDM.preTest(currentSignal != null);

        return currentSignal;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: "

            + "CurrentBloodPressure Level not Intialized");

    }

    try{

        VDM.invTest(this);

    } catch(Exception e) {

        System.err.println("Postcondition test failed: "
```

Ambulatory Blood Pressure Monitoring System

```
+ e.getMessage());  
  
}  
  
return null;  
  
}  
  
public BloodPressureLevel determineBloodPressureLevel() {  
  
    if (systolic < 90 && diastolic < 60) {  
  
        try {  
  
            VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);  
  
            currentBloodPressureLevel = BloodPressureLevel.LOW_BLOOD_PRESSURE;  
  
        } catch (VDM.VDMPreconditionException e) {  
  
            System.err.println("Precondition test failed: " + e.getMessage());  
  
        }  
  
        try {  
  
            VDM.invTest(this);  
  
        } catch (Exception e) {  
  
            System.err.println("Postcondition test failed: "  
  
                + e.getMessage());  
  
        }  
  
        return currentBloodPressureLevel;  
  
    } else if (systolic >= 90 && systolic <= 120 && diastolic >= 60  
  
        && diastolic <= 80) {  
  
        try {  
  
            VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);  
  
            currentBloodPressureLevel = BloodPressureLevel.NORMAL;  
  
        } catch (VDM.VDMPreconditionException e) {  
  
            System.err.println("Precondition test failed: " + e.getMessage());  
  
        }  
  
    }  
  
}
```

Ambulatory Blood Pressure Monitoring System

```
}

try{

VDM.invTest(this);

} catch(Exception e) {

    System.err.println("Postcondition test failed: "

        + e.getMessage());

}

return currentBloodPressureLevel;

} else if (systolic > 120 && systolic <= 130 && diastolic <= 80) {

try {

    VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);

    currentBloodPressureLevel = BloodPressureLevel.ELEVATED;

} catch (VDM.VDMPreconditionException e) {

    System.err.println("Precondition test failed: " + e.getMessage());

}

try{

VDM.invTest(this);

} catch(Exception e) {

    System.err.println("Postcondition test failed: "

        + e.getMessage());

}

return currentBloodPressureLevel;

} else if (systolic > 130 && systolic <= 140 && diastolic > 80

    && diastolic <= 90) {

try {

    VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);
```

Ambulatory Blood Pressure Monitoring System

```
        currentBloodPressureLevel = BloodPressureLevel.HBP_STAGE_1;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: " + e.getMessage());

    }

    try{

        VDM.invTest(this);

    } catch(Exception e) {

        System.err.println("Postcondition test failed: "

            + e.getMessage());

    }

    return currentBloodPressureLevel;

} else if (systolic > 140 && diastolic > 90) {

    try {

        VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);

        currentBloodPressureLevel = BloodPressureLevel.HYPERTENSIVE_CRISIS;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: " + e.getMessage());

    }

    try{

        VDM.invTest(this);

    } catch(Exception e) {

        System.err.println("Postcondition test failed: "

            + e.getMessage());

    }

    return currentBloodPressureLevel;

} else {
```

Ambulatory Blood Pressure Monitoring System

```
try {  
    VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);  
    currentBloodPressureLevel = BloodPressureLevel.HBP_STAGE_2;  
} catch (VDM.VDMPreconditionException e) {  
    System.err.println("Precondition test failed: " + e.getMessage());  
}  
  
try{  
    VDM.invTest(this);  
} catch(Exception e) {  
    System.err.println("Postcondition test failed: "  
        + e.getMessage());  
}  
  
return currentBloodPressureLevel;  
}  
  
}  
  
public SIGNAL determineSignal() {  
    if (systolic < 90 && diastolic < 60) {  
        try {  
            VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);  
            currentSignal = SIGNAL.VIBRATE_AND_ALERT;  
        } catch (VDM.VDMPreconditionException e) {  
            System.err.println("Precondition test failed: " + e.getMessage());  
        }  
        try{  
            VDM.invTest(this);  
        } catch(Exception e) {
```


Ambulatory Blood Pressure Monitoring System

```
System.err.println("Postcondition test failed: "
    + e.getMessage());
}

return currentSignal;
} else if (systolic >= 90 && systolic <= 120 && diastolic >= 60
    && diastolic <= 80) {
    try {
        VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);
        currentSignal = SIGNAL.SILENT;
    } catch (VDM.VDMPreconditionException e) {
        System.err.println("Precondition test failed: " + e.getMessage());
    }
    try{
        VDM.invTest(this);
    } catch(Exception e) {
        System.err.println("Postcondition test failed: "
            + e.getMessage());
    }
    return currentSignal;
} else if (systolic > 120 && systolic <= 130 && diastolic <= 80) {
    try {
        VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);
        currentSignal = SIGNAL.VIBRATE;
    } catch (VDM.VDMPreconditionException e) {
        System.err.println("Precondition test failed: " + e.getMessage());
    }
}
```

Ambulatory Blood Pressure Monitoring System

```
try{

    VDM.invTest(this);

} catch(Exception e) {

    System.err.println("Postcondition test failed: "

        + e.getMessage());

}

return currentSignal;

} else if (systolic > 130 && systolic <= 140 && diastolic > 80

    && diastolic <= 90) {

    try {

        VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);

        currentSignal = SIGNAL.VIBRATE_AND_ALERT;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: " + e.getMessage());

    }

    try{

        VDM.invTest(this);

    } catch(Exception e) {

        System.err.println("Postcondition test failed: "

            + e.getMessage());

    }

    return currentSignal;

} else if (systolic > 140 && diastolic > 90) {

    try {

        VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);

        currentSignal = SIGNAL.VIBRATE_AND_ALERT;
```

Ambulatory Blood Pressure Monitoring System

```
    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: " + e.getMessage());

    }

    try{

        VDM.invTest(this);

    } catch (Exception e) {

        System.err.println("Postcondition test failed: "

            + e.getMessage());

    }

    return currentSignal;

} else {

    try {

        VDM.preTest(inRange(systolic) && inRange(diastolic) && systolic!=0 && diastolic!=0);

        currentSignal = SIGNAL.VIBRATE_AND_ALERT;

    } catch (VDM.VDMPreconditionException e) {

        System.err.println("Precondition test failed: " + e.getMessage());

    }

    try{

        VDM.invTest(this);

    } catch (Exception e) {

        System.err.println("Postcondition test failed: "

            + e.getMessage());

    }

    return currentSignal;

}

}
```

Ambulatory Blood Pressure Monitoring System

```
public int getSystolic() {  
  
    try {  
  
        VDM.preTest(systolic!=0);  
  
    } catch (VDM.VDMPreconditionException e) {  
  
        System.err.println("Precondition test failed: " + "Stored Systolic Measure is out of Bound Or not initialized");  
  
    }  
  
    try{  
  
        VDM.invTest(this);  
  
    } catch(Exception e) {  
  
        System.err.println("Postcondition test failed: "  
            + e.getMessage());  
  
    }  
  
    return systolic;  
  
}  
  
public int getDiastolic() {  
  
    try {  
  
        VDM.preTest(diastolic!=0);  
  
    } catch (VDM.VDMPreconditionException e) {  
  
        System.err.println("Precondition test failed: " + "Stored Diastolic Measure is out of Bound Or not initialized");  
  
    }  
  
    try{  
  
        VDM.invTest(this);  
  
    } catch(Exception e) {  
  
        System.err.println("Postcondition test failed: "  
            + e.getMessage());  
  
    }  
  
}
```

Ambulatory Blood Pressure Monitoring System

```
    return diastolic;  
}  
  
}
```

VI. TEST CASES SUMMARY

Test Case ID	Test Case	Input	Expected Output	Actual Output
01	InitializingSystolicPressure_InvalidInput	pressureIn = 500	Exception Thrown: Pre-condition failed. And systolic=0	Exception: Pre-condition failed. systolic=0
02	InitializingSystolicPressure_ValidInput	pressureIn=150	systolic=150	systolic=150
03	InitializingDiastolicPressure_InvalidInput	pressureIn=450	Exception Thrown: Pre-condition failed. And diastolic=0	Exception: Pre-condition failed. diastolic=0
04	InitializingDiastolicPressure_ValidInput	pressureIn=100	diastolic=100	diastolic=100
05	ChangingSystolicPressure_InvalidInput	pressureIn=410	Exception Thrown: Pre-condition Failed. And Systolic does not change.	Exception: Pre-condition failed.
06	ChangingSystolicPressure_ValidInput	pressureIn=200	systolic=200	systolic=200
07	ChangingDiastolic	pressureIn=460	Exception	Exception:

Ambulatory Blood Pressure Monitoring System

	Pressure_InvalidInput		Thrown: Pre-condition Failed. And Diastolic does not change.	Pre-condition failed.
08	ChangingDiastolic Pressure_ValidInput	pressureIn=160	diastolic=160	diastolic=160
09	GetBloodPressure Level_AfterInitialization	Systolic and diastolic are initialized. And Systolic=100 Diastolic=70	BloodPressure=<NORMAL>	BloodPressure=<NORMAL>
10	GetBloodPressure Level_BeforeInitialization	Systolic and diastolic are not initialized.	Exception Thrown: Pre-condition Failed.	Exception: Pre-condition failed.
11	GetSystolicPressure_AfterInitialization	Systolic is initialized. And systolic=70	pressureOut=70	pressureOut=70
12	GetSystolicPressure_BeforeInitialization	Systolic is not initialized.	Exception Thrown: Pre-condition Failed.	Exception: Pre-condition failed.
13	GetDiastolicPressure_AfterInitialization	Diastolic is initialized. And diastolic=70	pressureOut=70	pressureOut=70
14	GetDiastolicPressure_BeforeInitialization	Diastolic is not initialized.	Exception Thrown: Pre-condition Failed.	Exception: Pre-condition failed.
15	DetermineBloodPressure_LOWBP_AND_DetermineSig	systolic=85 diastolic=55	bloodPressureOutput = <LOW_BLOOD	bloodPressureOutput = LOW_BLOOD_

Ambulatory Blood Pressure Monitoring System

	nal_<ALERT_AND_VIBRATE>		_PRESSURE>	PRESSURE
16	DetermineBloodPressure_NORMAL_AND_DetermineSignal_<SILENT>	systolic=100 diastolic=70	bloodPressureOutput = <NORMAL> AND signalOut=<SILENT>	bloodPressureOutput = NORMAL signalOut=<SILENT>
17	DetermineBloodPressure_ELEVATED_AND_DetermineSignal_<VIBRATE>	systolic=121 diastolic=78	bloodPressureOutput = <ELEVATED> AND signalOut=<VIBRATE>	bloodPressureOutput = ELEVATED signalOut=<VIBRATE>
18	DetermineBloodPressure_HBPSTAGE1_AND_DetermineSignal_<ALERT_AND_VIBRATE>	systolic=135 diastolic=85	bloodPressureOutput = <HBP_STAGE1> AND signalOut=<ALERT_AND_VIBRATE>	bloodPressureOutput = HBP_STAGE1 signalOut=<ALERT_AND_VIBRATE>
19	DetermineBloodPressure_HBPSTAGE2_AND_DetermineSignal_<ALERT_AND_VIBRATE>	systolic=145 diastolic=85	bloodPressureOutput = <HBP_STAGE2> AND signalOut=<ALERT_AND_VIBRATE>	bloodPressureOutput = HBP_STAGE2 signalOut=<ALERT_AND_VIBRATE>
20	DetermineBloodPressure_HYPERTENSIVECRISIS_AND_DetermineSignal_<ALERT_AND_VIBRATE>	systolic=145 diastolic=100	bloodPressureOutput = <HYPERTENSIVE_CRISIS> And signalOut=<ALERT_AND_VIBRATE>	bloodPressureOutput = HYPERTENSIVE_CRISIS signalOut=<ALERT_AND_VIBRATE>

Ambulatory Blood Pressure Monitoring System

VII. TESTER CLASS JAVA CODE :

```
import java.util.Scanner;

import ABPM.AmbulatoryBloodPressureController;

public class Tester {
    public static void main(String[] args) {
        Scanner scanner = new Scanner(System.in);

        AmbulatoryBloodPressureController abpc =
            new AmbulatoryBloodPressureController("John", "Doe", 62);

        int choice;

        do {
            System.out.println("\n----- Tester Menu -----");
            System.out.println("1. InitializingSystolicPressure_InvalidInput");
            System.out.println("2. InitializingSystolicPressure_ValidInput");
            System.out.println("3. InitializingDiastolicPressure_InvalidInput");
            System.out.println("4. InitializingDiastolicPressure_ValidInput");
            System.out.println("5. ChangingSystolicPressure_InvalidInput");
            System.out.println("6. ChangingSystolicPressure_ValidInput");
            System.out.println("7. ChangingDiastolicPressure_InvalidInput");
            System.out.println("8. ChangingDiastolicPressure_ValidInput");
            System.out.println("9. GetBloodPressureLevel_BeforeInitialization");
            System.out.println("10. GetSystolicPressure_AfterInitilization");
            System.out.println("11. GetSystolicPressure_BeforeInitialization");
            System.out.println("12. GetDiastolicPressure_AfterInitilization");
            System.out.println("13. GetDiastolicPressure_BeforeInitialization");
            System.out.println("14. DetermineBloodPressure_LOWBP &&
DetermineSignal_<ALERT_AND_VIBRATE>");
            System.out.println("15. DetermineBloodPressure_NORMAL &&
DetermineSignal_<SILENT>");
            System.out.println("16. DetermineBloodPressure_ELEVATED &&
DetermineSignal_<VIBRATE>");
            System.out.println("17. DetermineBloodPressure_HBPSTAGE1 &&
DetermineSignal_<ALERT_AND_VIBRATE>");
            System.out.println("18. DetermineBloodPressure_HBPSTAGE2 &&
DetermineSignal_<ALERT_AND_VIBRATE>");
```


Ambulatory Blood Pressure Monitoring System

```
System.out.println("19. DetermineBloodPressure_HYPERTENSIVECRISIS &&
DetermineSignal_<ALERT_AND_VIBRATE>");
System.out.println("0. Exit");

System.out.print("Enter your choice: ");
choice = scanner.nextInt();

switch (choice) {
    case 1:
        System.out.println("Expected Output: Exception thrown and Systolic set to 0");
        abpc.setInitialSystolicPressure(500);
        System.out.println("Actual Output: "+abpc.getSystolic());
        break;

    case 2:
        System.out.println("Expected Output: Systolic set to 150");
        abpc.setInitialSystolicPressure(150);
        System.out.println("Actual Output: "+abpc.getSystolic());
        break;

    case 3:
        System.out.println("Expected Output: Exception thrown and Diastolic set to 0");
        abpc.setInitialDiastolicPressure(450);
        System.out.println("Actual Output: "+abpc.getDiastolic());
        break;

    case 4:
        System.out.println("Expected Output: Diastolic set to 100");
        abpc.setInitialDiastolicPressure(100);
        System.out.println("Actual Output: "+abpc.getDiastolic());
        break;

    case 5:
        System.out.println("Expected Output: Exception Thrown: Pre-condition Failed." +
            "Systolic does not Changed (0)");
        abpc.changeSystolicPressure(410);
        System.out.println("Actual Output: "+abpc.getSystolic());
        break;

    case 6:
```

Ambulatory Blood Pressure Monitoring System

```
System.out.println("Expected Output: Systolic = 200");
abpc.changeSystolicPressure(200);
System.out.println("Actual Output: "+abpc.getSystolic());
break;
```

case 7:

```
System.out.println("Expected Output: Exception Thrown: Pre-condition Failed." +
    "Diastolic does not Changed (0)");
abpc.changeDiastolicPressure(460);
System.out.println("Actual Output: "+abpc.getDiastolic());
break;
```

case 8:

```
System.out.println("Expected Output: Systolic = 160");
abpc.changeDiastolicPressure(160);
System.out.println("Actual Output: "+abpc.getDiastolic());
break;
```

case 9:

```
                AmbulatoryBloodPressureController newone =new
AmbulatoryBloodPressureController("John", "Doe", 62);
System.out.println("Expected Output: Exception Thrown: Pre-condition Failed.");
System.out.println("Actual Output: ");
newone.getCurrentBloodPressureLevel();
break;
```

case 10:

```
                System.out.println("Expected Output: Systolic is initialized." + "And" +
"systolic=70" + "");
abpc.setInitialSystolicPressure(70);
System.out.println("Actual Output: "+abpc.getSystolic());
break;
```

case 11:

```
                AmbulatoryBloodPressureController newtwo =new
AmbulatoryBloodPressureController("John", "Doe", 62);
System.out.println("Exception Thrown: Pre-condition Failed.");
System.out.println("Actual Output: "+newtwo.getSystolic());
break;
```

Ambulatory Blood Pressure Monitoring System

```
case 12:
    System.out.println("Expected Output: Diastolic is initialized." + "And" + "Diastolic
=70" + "");
    abpc.setInitialDiastolicPressure(70);
    System.out.println("Actual Output: "+abpc.getDiastolic());
    break;

case 13:
    AmbulatoryBloodPressureController newthree =new
AmbulatoryBloodPressureController("John", "Doe", 62);
    System.out.println("Expected Output: Exception Thrown: Pre-condition Failed.");
    System.out.println("Actual Output: "+newthree.getDiastolic());
    break;

case 14:
    System.out.println("Expected Output: bloodPressureOut =
<LOW_BLOOD_PRESSURE> and Signal: <ALERT_AND_VIBRATE>");
    abpc.setInitialSystolicPressure(150);
    abpc.setInitialDiastolicPressure(100);
    abpc.changeDiastolicPressure(55);
    abpc.changeSystolicPressure(85);
    System.out.println("Actual Output: "+abpc.determineBloodPressureLevel()+"
"+abpc.determineSignal());
    break;

case 15:
    System.out.println("Expected Output: bloodPressureOut = <NORMAL> and Signal:
<SILENT>");
    abpc.setInitialSystolicPressure(150);
    abpc.setInitialDiastolicPressure(100);
    abpc.changeDiastolicPressure(70);
    abpc.changeSystolicPressure(100);
    System.out.println("Actual Output: "+abpc.determineBloodPressureLevel()+"
"+abpc.determineSignal());
    break;

case 16:
    System.out.println("Expected Output: bloodPressureOut = <ELEVATED> and
Signal: <VIBRATE>");
    abpc.setInitialSystolicPressure(150);
```

Ambulatory Blood Pressure Monitoring System

```
        abpc.setInitialDiastolicPressure(100);
        abpc.changeDiastolicPressure(78);
        abpc.changeSystolicPressure(121);
        System.out.println("Actual Output: "+abpc.determineBloodPressureLevel()+"
"+abpc.determineSignal());
        break;

    case 17:
        System.out.println("Expected Output: bloodPressureOut = <HBP_STAGE1> and
Signal: <ALERT_AND_VIBRATE>");
        abpc.setInitialSystolicPressure(150);
        abpc.setInitialDiastolicPressure(100);
        abpc.changeDiastolicPressure(85);
        abpc.changeSystolicPressure(135);
        System.out.println("Actual Output: "+abpc.determineBloodPressureLevel()+"
"+abpc.determineSignal());
        break;

    case 18:
        System.out.println("Expected Output: bloodPressureOut = <HBP_STAGE2> and
Signal: <ALERT_AND_VIBRATE>");
        abpc.setInitialSystolicPressure(150);
        abpc.setInitialDiastolicPressure(100);
        abpc.changeDiastolicPressure(85);
        abpc.changeSystolicPressure(145);
        System.out.println("Actual Output: "+abpc.determineBloodPressureLevel()+"
"+abpc.determineSignal());
        break;

    case 19:
        System.out.println("Expected Output: bloodPressureOut =
<HYPERTENSIVECRISIS> and Signal: <ALERT_AND_VIBRATE>");
        abpc.setInitialSystolicPressure(150);
        abpc.setInitialDiastolicPressure(100);
        abpc.changeDiastolicPressure(100);
        abpc.changeSystolicPressure(145);
        System.out.println("Actual Output: "+abpc.determineBloodPressureLevel()+"
"+abpc.determineSignal());
        break;
```

Ambulatory Blood Pressure Monitoring System

```
        case 0:
            System.out.println("Shutting Down the device. Goodbye!");
            break;

        default:
            System.out.println("Invalid choice. Please try again.");
            break;
    }

    } while (choice != 0);
}
```